



TRPM1 Antibody, Biotin conjugated

Product Code	CSB-PA773791LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z4N2
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily M member 1 protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily M member 1 (Long transient receptor potential channel 1) (LTrpC1) (Melastatin-1), TRPM1, LTRPC1 MLSN MLSN1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRPM1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.