



TMEM150A Antibody, FITC conjugated

Product Code	CSB-PA774796LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86TG1
Immunogen	Recombinant Human Transmembrane protein 150A protein (24-75AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transmembrane protein 150A (Transmembrane protein 150), TMEM150A, TMEM150
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TMEM150A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.