



TIPARP Antibody, Biotin conjugated

Product Code	CSB-PA801800LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Z3E1
Immunogen	Recombinant Human TCDD-inducible poly [ADP-ribose] polymerase protein (74-178AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TCDD-inducible poly [ADP-ribose] polymerase (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 14) (ARTD14) (Poly [ADP-ribose] polymerase 7) (PARP-7), TIPARP, PARP7
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TIPARP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.