



MAML2 Antibody, Biotin conjugated

Product Code	CSB-PA809028LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZL2
Immunogen	Recombinant Human Mastermind-like protein 2 protein (347-506AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mastermind-like protein 2 (Mam-2), MAML2, KIAA1819
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MAML2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.