



ZDHHC14 Antibody, HRP conjugated

Product Code	CSB-PA811637LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZN3
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC14 protein (328-453AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable palmitoyltransferase ZDHHC14 (EC 2.3.1.225) (NEW1 domain-containing protein) (NEW1CP) (Zinc finger DHHC domain-containing protein 14) (DHHC-14), ZDHHC14
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZDHHC14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.