



KIRREL3 Antibody, Biotin conjugated

Product Code	CSB-PA811640LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IZU9
Immunogen	Recombinant Human Kin of IRRE-like protein 3 protein (200-343AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kin of IRRE-like protein 3 (Kin of irregular chiasm-like protein 3) (Nephrin-like protein 2) [Cleaved into: Processed kin of IRRE-like protein 3], KIRREL3, KIAA1867 NEPH2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KIRREL3