




PDILT Antibody

Product Code	CSB-PA818726LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N807
Immunogen	Recombinant Human Protein disulfide-isomerase-like protein of the testis protein (494-584AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase-like protein of the testis, PDILT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDILT

Image



Western Blot

Positive WB detected in: A549 whole cell lysate, Colo320 whole cell lysate

All lanes: PDILT antibody at 4.63µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 67 kDa

Observed band size: 67 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.