



PDILT Antibody, Biotin conjugated

Product Code	CSB-PA818726LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N807
Immunogen	Recombinant Human Protein disulfide-isomerase-like protein of the testis protein (494-584AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase-like protein of the testis, PDILT
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDILT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.