



FANCB Antibody, HRP conjugated

Product Code	CSB-PA818751LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NB91
Immunogen	Recombinant Human Fanconi anemia group B protein (431-574AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fanconi anemia group B protein (Protein FACB) (Fanconi anemia-associated polypeptide of 95 kDa) (FAAP95), FANCB
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FANCB