



FKBP10 Antibody, FITC conjugated

Product Code	CSB-PA822180LC01HU
Abbreviation	Peptidyl-prolyl cis-trans isomerase FKBP10
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96AY3
Immunogen	Recombinant Human Peptidyl-prolyl cis-trans isomerase FKBP10 protein (459-579AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	PPIases accelerate the folding of proteins during protein synthesis.
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-prolyl cis-trans isomerase FKBP10 (PPIase FKBP10) (EC 5.2.1.8) (65 kDa FK506-binding protein) (65 kDa FKBP) (FKBP-65) (FK506-binding protein 10) (FKBP-10) (Immunophilin FKBP65) (Rotamase), FKBP10, FKBP65
Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FKBP10