



SPPL2B Antibody, Biotin conjugated

Product Code	CSB-PA823457LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCT7
Immunogen	Recombinant Human Signal peptide peptidase-like 2B protein (492-592AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Signal peptide peptidase-like 2B (SPP-like 2B) (SPPL2b) (EC 3.4.23.-) (Intramembrane protease 4) (IMP-4) (Presenilin homologous protein 4) (PSH4) (Presenilin-like protein 1), SPPL2B, IMP4 KIAA1532 PSL1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	SPPL2B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.