



LMO7 Antibody, HRP conjugated

Product Code	CSB-PA823897LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WWI1
Immunogen	Recombinant Human LIM domain only protein 7 protein (1301-1465AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LIM domain only protein 7 (LMO-7) (F-box only protein 20) (LOMP), LMO7, FBX20 FBXO20 KIAA0858
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LMO7