



CTAGE1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA836295LD01HU |
| Abbreviation | cTAGE family member 2 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96RT6 |
| Immunogen | Recombinant Human cTAGE family member 2 protein (471-595AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | cTAGE family member 2 (Protein cTAGE-2) (Cancer/testis antigen 21.2) (CT21.2), CTAGE1, CTAGE2 |
| Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | CTAGE1 |