



SSH1 Antibody, HRP conjugated

Product Code	CSB-PA837881LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WYL5
Immunogen	Recombinant Human Protein phosphatase Slingshot homolog 1 protein (928-1049AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein phosphatase Slingshot homolog 1 (EC 3.1.3.16) (EC 3.1.3.48) (SSH-like protein 1) (SSH-1L) (hSSH-1L), SSH1, KIAA1298 SSH1L
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SSH1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.