



SNED1 Antibody, HRP conjugated

Product Code	CSB-PA840582LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TER0
Immunogen	Recombinant Human Sushi, nidogen and EGF-like domain-containing protein 1 protein (1200-1341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sushi, nidogen and EGF-like domain-containing protein 1 (Insulin-responsive sequence DNA-binding protein 1) (IRE-BP1), SNED1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SNED1