




POF1B Antibody

Product Code	CSB-PA840988LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVV4
Immunogen	Recombinant Human Protein POF1B protein (532-589AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein POF1B (Premature ovarian failure protein 1B), POF1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POF1B

Image



Western Blot
 Positive WB detected in: K562 whole cell lysate,
 Rat lung tissue
 All lanes: POF1B antibody at 3.4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 69, 35 kDa
 Observed band size: 69 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.