



POF1B Antibody, FITC conjugated

Product Code	CSB-PA840988LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVV4
Immunogen	Recombinant Human Protein POF1B protein (532-589AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein POF1B (Premature ovarian failure protein 1B), POF1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POF1B