



POF1B Antibody, Biotin conjugated

Product Code	CSB-PA840988LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WVV4
Immunogen	Recombinant Human Protein POF1B protein (532-589AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein POF1B (Premature ovarian failure protein 1B), POF1B
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	POF1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.