



MARVELD3 Antibody, HRP conjugated

Product Code	CSB-PA842633LB01HU
Abbreviation	MARVEL domain-containing protein 3
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96A59
Immunogen	Recombinant Human MARVEL domain-containing protein 3 protein (1-189AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Relevance	As a component of tight junctions, plays a role in paracellular ion conductivity.
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MARVEL domain-containing protein 3, MARVELD3, MRVLDC3
Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MARVELD3