



ZFR Antibody, Biotin conjugated

Product Code	CSB-PA842700LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96KR1
Immunogen	Recombinant Human Zinc finger RNA-binding protein (150-268AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger RNA-binding protein (hZFR) (M-phase phosphoprotein homolog), ZFR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZFR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.