



IFT140 Antibody, HRP conjugated

Product Code	CSB-PA842756LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RY7
Immunogen	Recombinant Human Intraflagellar transport protein 140 homolog protein (1291-1462AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 140 homolog (WD and tetratricopeptide repeats protein 2), IFT140, KIAA0590 WDTC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	IFT140
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.