



POC5 Antibody, Biotin conjugated

Product Code	CSB-PA843309LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NA72
Immunogen	Recombinant Human Centrosomal protein POC5 protein (373-485AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein POC5 (Protein of centriole 5) (hPOC5), POC5, C5orf37
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	POC5