



NKX2-3 Antibody, Biotin conjugated

Product Code	CSB-PA848391LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TAU0
Immunogen	Recombinant Human Homeobox protein Nkx-2.3 protein (52-148AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Nkx-2.3 (Homeobox protein NK-2 homolog C), NKX2-3, NKX23 NKX2C
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NKX2-3