



MFSD4A Antibody, FITC conjugated

Product Code	CSB-PA850797LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N468
Immunogen	Recombinant Human Major facilitator superfamily domain-containing protein 4A protein (160-220AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major facilitator superfamily domain-containing protein 4A (Major facilitator superfamily domain-containing protein 4), MFSD4A, MFSD4
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MFSD4A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.