



SNX33 Antibody, Biotin conjugated

Product Code	CSB-PA851969LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WV41
Immunogen	Recombinant Human Sorting nexin-33 protein (348-516AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sorting nexin-33 (SH3 and PX domain-containing protein 3), SNX33, SH3PX3 SH3PXD3C SNX30
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SNX33
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.