



TOP1MT Antibody, FITC conjugated

Product Code	CSB-PA853372LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969P6
Immunogen	Recombinant Human DNA topoisomerase I, mitochondrial protein (270-400AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA topoisomerase I, mitochondrial (TOP1mt) (EC 5.99.1.2), TOP1MT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	TOP1MT