



SMC1B Antibody, Biotin conjugated

Product Code	CSB-PA854104LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NDV3
Immunogen	Recombinant Human Structural maintenance of chromosomes protein 1B protein (954-1060AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Structural maintenance of chromosomes protein 1B (SMC protein 1B) (SMC-1-beta) (SMC-1B), SMC1B, SMC1L2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SMC1B