



MAGI1 Antibody, FITC conjugated

Product Code	CSB-PA857014LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96QZ7
Immunogen	Recombinant Human Membrane-associated guanylate kinase, WW and PDZ domain-containing protein 1 protein (1234-1367AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Membrane-associated guanylate kinase, WW and PDZ domain-containing protein 1 (Atrophin-1-interacting protein 3) (AIP-3) (BAI1-associated protein 1) (BAP-1) (Membrane-associated guanylate kinase inverted 1) (MAGI-1) (Trinucleotide repeat-containing gene 19 protein) (WW domain-containing protein 3) (WWP3), MAGI1, AIP3 BAIAP1 BAP1 TNRC19
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MAGI1