



FEV Antibody, FITC conjugated

Product Code	CSB-PA858717LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99581
Immunogen	Recombinant Human Protein FEV protein (1-97AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FEV (Fifth Ewing variant protein) (PC12 ETS domain-containing transcription factor 1) (PC12 ETS factor 1) (Pet-1), FEV, PET1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	FEV
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.