



FSD1 Antibody, HRP conjugated

Product Code	CSB-PA861136LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV5
Immunogen	Recombinant Human Fibronectin type III and SPRY domain-containing protein 1 protein (251-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin type III and SPRY domain-containing protein 1 (MID1-related protein 1) (Microtubule-associated protein GLFND), FSD1, GLFND MIR1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FSD1