



FSD1 Antibody, HRP conjugated

Product Code	CSB-PA861136LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTV5
Immunogen	Recombinant Human Fibronectin type III and SPRY domain-containing protein 1 protein (251-370AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibronectin type III and SPRY domain-containing protein 1 (MID1-related protein 1) (Microtubule-associated protein GLFND), FSD1, GLFND MIR1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FSD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.