



NLN Antibody, HRP conjugated

Product Code	CSB-PA863161LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYT8
Immunogen	Recombinant Human Neurolysin, mitochondrial protein (109-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Neurolysin, mitochondrial (EC 3.4.24.16) (Angiotensin-binding protein) (Microsomal endopeptidase) (MEP) (Mitochondrial oligopeptidase M) (Neurotensin endopeptidase), NLN, AGTBP KIAA1226
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NLN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.