




PCGF6 Antibody

Product Code	CSB-PA871612LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYE7
Immunogen	Recombinant Human Polycomb group RING finger protein 6 protein (9-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb group RING finger protein 6 (Mel18 and Bmi1-like RING finger) (RING finger protein 134), PCGF6, MBLR RNF134
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PCGF6

Image



Western Blot

Positive WB detected in: PC-3 whole cell lysate, MCF-7 whole cell lysate

All lanes: PCGF6 antibody at 3.1µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 40, 31 kDa

Observed band size: 40 kDa