



PCGF6 Antibody, Biotin conjugated

Product Code	CSB-PA871612LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYE7
Immunogen	Recombinant Human Polycomb group RING finger protein 6 protein (9-70AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb group RING finger protein 6 (Mel18 and Bmi1-like RING finger) (RING finger protein 134), PCGF6, MBLR RNF134
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	PCGF6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.