



INPP5E Antibody, HRP conjugated

Product Code	CSB-PA873641LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NRR6
Immunogen	Recombinant Human 72 kDa inositol polyphosphate 5-phosphatase protein (1-133AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	72 kDa inositol polyphosphate 5-phosphatase (EC 3.1.3.36) (Phosphatidylinositol 4,5-bisphosphate 5-phosphatase) (Phosphatidylinositol polyphosphate 5-phosphatase type IV), INPP5E
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5E
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.