



SLC4A5 Antibody, Biotin conjugated

Product Code	CSB-PA874837LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BY07
Immunogen	Recombinant Human Electrogenic sodium bicarbonate cotransporter 4 protein (1071-1137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Electrogenic sodium bicarbonate cotransporter 4 (NBCe2) (Solute carrier family 4 member 5), SLC4A5, NBC4
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SLC4A5