



PEAK1 Antibody, FITC conjugated

Product Code	CSB-PA875688LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H792
Immunogen	Recombinant Human Pseudopodium-enriched atypical kinase 1 protein (715-881AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pseudopodium-enriched atypical kinase 1 (EC 2.7.10.2) (Sugen kinase 269) (Tyrosine-protein kinase SgK269), PEAK1, KIAA2002
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PEAK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.