



DOHH Antibody, FITC conjugated

Product Code	CSB-PA880104LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BU89
Immunogen	Recombinant Human Deoxyhypusine hydroxylase protein (79-159AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Deoxyhypusine hydroxylase (hDOHH) (EC 1.14.99.29) (Deoxyhypusine dioxygenase) (Deoxyhypusine monooxygenase) (HEAT-like repeat-containing protein 1), DOHH, HLRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	DOHH
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.