



TBCD Antibody, Biotin conjugated

Product Code	CSB-PA883599LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BTW9
Immunogen	Recombinant Human Tubulin-specific chaperone D protein (2-96AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin-specific chaperone D (Beta-tubulin cofactor D) (tfcD) (SSD-1) (Tubulin-folding cofactor D), TBCD, KIAA0988 SSD1 TFCD
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TBCD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.