



# BACH2 Antibody, Biotin conjugated

Product Code	CSB-PA883639LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYV9
Immunogen	Recombinant Human Transcription regulator protein BACH2 protein (136-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription regulator protein BACH2 (BTB and CNC homolog 2), BACH2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	BACH2