



BACH2 Antibody, Biotin conjugated

Product Code	CSB-PA883639LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYV9
Immunogen	Recombinant Human Transcription regulator protein BACH2 protein (136-235AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription regulator protein BACH2 (BTB and CNC homolog 2), BACH2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	BACH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.