



TRPM3 Antibody, Biotin conjugated

Product Code	CSB-PA884520LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9HCF6
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily M member 3 protein (1606-1713AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily M member 3 (Long transient receptor potential channel 3) (LTrpC-3) (LTrpC3) (Melastatin-2) (MLSN2), TRPM3, KIAA1616 LTRPC3
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TRPM3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.