



FAM83D Antibody, FITC conjugated

Product Code	CSB-PA887979LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4H8
Immunogen	Recombinant Human Protein FAM83D protein (339-470AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM83D (Spindle protein CHICA), FAM83D, C20orf129
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	FAM83D
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.