



KCNK9 Antibody, FITC conjugated

Product Code	CSB-PA889068LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPC2
Immunogen	Recombinant Human Potassium channel subfamily K member 9 protein (250-362AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Potassium channel subfamily K member 9 (Acid-sensitive potassium channel protein TASK-3) (TWIK-related acid-sensitive K(+) channel 3) (Two pore potassium channel KT3.2) (Two pore K(+) channel KT3.2), KCNK9, TASK3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNK9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.