



POLE4 Antibody, HRP conjugated

Product Code	CSB-PA889087LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR33
Immunogen	Recombinant Human DNA polymerase epsilon subunit 4 protein (8-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	DNA polymerase epsilon subunit 4 (EC 2.7.7.7) (DNA polymerase II subunit 4) (DNA polymerase epsilon subunit p12), POLE4
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	POLE4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.