



INTU Antibody, Biotin conjugated

Product Code	CSB-PA891544LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULD6
Immunogen	Recombinant Human Protein inturned protein (6-139AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein inturned (Inturned planar cell polarity effector homolog) (PDZ domain-containing protein 6), INTU, KIAA1284 PDZD6 PDZK6
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	INTU
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.