



INO80 Antibody, Biotin conjugated

Product Code	CSB-PA892466LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ULG1
Immunogen	Recombinant Human Chromatin-remodeling ATPase INO80 protein (1365-1495AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chromatin-remodeling ATPase INO80 (hINO80) (DNA helicase-related INO80 complex homolog 1) (DNA helicase-related protein INO80) (INO80 complex subunit A), INO80, INO80A, INOC1, KIAA1259
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	INO80
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.