



TNN Antibody, Biotin conjugated

Product Code	CSB-PA892498LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UQP3
Immunogen	Recombinant Human Tenascin-N protein (441-556AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tenascin-N (TN-N) (Tenascin-W) (TN-W), TNN, TNW
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TNN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.