



NUBP2 Antibody, Biotin conjugated

Product Code	CSB-PA896918LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Y2
Immunogen	Recombinant Human Cytosolic Fe-S cluster assembly factor NUBP2 protein (54-170AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic Fe-S cluster assembly factor NUBP2 (Nucleotide-binding protein 2) (NBP 2), NUBP2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	NUBP2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.