



RAPGEF2 Antibody, HRP conjugated

Product Code	CSB-PA897493LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y4G8
Immunogen	Recombinant Human Rap guanine nucleotide exchange factor 2 protein (1393-1498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rap guanine nucleotide exchange factor 2 (Cyclic nucleotide ras GEF) (CNrasGEF) (Neural RAP guanine nucleotide exchange protein) (nRap GEP) (PDZ domain-containing guanine nucleotide exchange factor 1) (PDZ-GEF1) (RA-GEF-1) (Ras/Rap1-associating GEF-1), RAPGEF2, KIAA0313 NRAPGEP PDZGEF1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RAPGEF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.